

## MAGNETIC FIELD IN A SOLENOID

### 1. EXPERIMENT OBJECTIVES

- 1.1 Study the magnitude of the magnetic field strength inside and outside the solenoid .
- 1.2 Determine the relationship between magnetic field strength and electric current,
- 1.3 Determine the relationship between magnetic field strength and coil density,
- 1.4 Determine constant permeability magnetic  $\mu_0$  .

### 2. TOOLS USED



Gambar 5.1. Tools used.

- 2.1 LabQuest2 interface device,
- 2.2 DC- variable power supply ,
- 2.3 Rheostat ,
- 2.4 Amperemeter,

- 2.5 Magnetic field sensor,
- 2.6 Solenoid standard board,
- 2.7 Solenoid and board,
- 2.8 Connecting wires.

### 3. BASIC CONCEPTS

A solenoid can be made by winding a wire around a tube in a certain number of turns. If an electric current flows through the solenoid, a magnetic field will be generated inside the solenoid. The magnitude of the magnetic field is expressed by the equation:

$$\mathbf{B} = \mu_0 n i \quad (5.1)$$

Where:

- B : magnetic field strength (Gauss / Tesla)
- $\mu_0$  : permeability constant (Tesla-meter/Ampere)
- n : number of turns per unit length of solenoid (1/m)
- i : electric current (Ampere)

## 4. EXPERIMENTAL STEPS

### 4.1 Experiment Tools Set Up

- 4.1.1 Set magnetic field sensor to *LabQuest* interface on Channel 1. set up sensor sensitivity to position *High*.
- 4.1.2 Stretch the Slinky solenoid to 1 m in length. Arrange so that the distance between the coils is approximately the same length.
- 4.1.3 Make a series circuit of the power supply, ammeter and rheostat. Connect the circuit to the solenoid.
- 4.1.4 Set the length of the solenoid to 1 meter. Place the sensor in the middle of the solenoid.
- 4.1.5 Turn on supply Power with voltage 10 V and set current Suite up to 500 mA.
- 4.1.6 Click the *Zero button* to ignore the influence of environmental magnetic fields.
- 4.1.7 Observe the magnitude of the magnetic field in *LabQuest*. Rotate the sensor so that the part marked white faces parallel to the solenoid axis. This is the most sensitive sensor position.
- 4.1.8 Press the button *play* to retrieve data. If the data shown by the Tesla Meter is stable enough, press the *Stop button* and record the value .
- 4.1.9 Reverse the sensor position, then record the measurement value.

### 4.2 Observe Strong Magnetic Field Relationship in Solenoid with Electric current

- 4.2.1 Set the length of the solenoid to 1 meter. Place the sensor in the middle of the solenoid.
- 4.2.2 Click the *Zero button* to ignore the influence of environmental magnetic fields.
- 4.2.3 Turn on supply power, set current circuit at 100 mA.
- 4.2.4 Press the button *play* to retrieve data. After data recording stops, record the average of measured value (using function *Analyze a Statistics*)
- 4.2.5 Repeat experimental steps 3 and 4 above for each 100 mA increase in current until a current of 700 mA is reached. (*consult with a practicum assistant*).
- 4.2.6 Repeat the experimental steps with measurements starting from a current of 700 mA to 100 mA (difference of 100 mA).
- 4.2.7 Tabulate the data that has been obtained in Table 5.1.

**Table 5.1.** Data from magnetic field measurements of the current applied to the solenoid.

length:				
No	Measurement goes up		Measurement down	
	$I(A)$	$B (mT)$	$I(A)$	$B (mT)$
1	0.1		0.7	
2	0.2		0.6	
3	0.3		0.5	
...	...		...	
7	0.7		0.1	

### 4.3 Observe Strong Magnetic Field Relationship in Solenoid with Coil Density

- 4.3.1 Set the solenoid length to 1.6 meters. Place the sensor in the middle of the solenoid.
- 4.3.2 Click the *Zero button* to ignore the influence of environmental magnetic fields.


- 4.3.3 Turn on supply power , set current Suite up to 700 mA,
- 4.3.4 Press the button *play* for retrieve data. After data recording stops, record the average of measured value ( using function *Analyze a Statistics* )
- 4.3.5 Repeat steps 1 to 4 for different numbers of coils per meter at several points. (consult with practicum assistant ).
- 4.3.6 Tabulate the data that has been obtained in Table 5.2.

**Table 5.2.** Data from magnetic field measurements on the number of coils per unit length of the solenoid.

Solenoid length:				
No	Measurement docked		Measurement stretch	
	$N/l$ ( turns / m )	$B$ ( mT )	$N/l$ ( turns / m )	$B$ ( mT )
1				
2				
3				
...				

#### 4.4 Observe the distribution of magnetic field values inside and outside the solenoid

- 4.4.1 Adjust the length of the solenoid until it reaches 1 meter.
- 4.4.2 Take 9 samples of magnetic field measurement points inside the solenoid as in Figure 5.2 without clicking the Zero button first .
- 4.4.3 Measure the magnetic field outside the solenoid, namely on each side of the solenoid .
- 4.4.4 Apply a current of 700 mA to the solenoid.
- 4.4.5 Record the measurement values according to the specified points.



**Figure 5.2** Location of measurement.

## 5. ANALYSIS AND DISCUSSIONS

- 5.1 Experiment with the tool setup section.
  - 5.1.1 Explain why it is necessary to use the zero button!
  - 5.1.2 Explain why the measured value of B changes sign (+, -) when it is reversed!

- 5.2 Experiment of magnetic field strength in the solenoid with current electricity relationship.
  - 5.2.1 Create a graph between the electric current as the Y axis and the magnetic field strength as the X axis (first convert the units to SI form)
  - 5.2.2 How is the magnetic field in solenoid a related to electric current ?
  - 5.2.3 Perform linear regression, then determine  $\mu_0$ .
  - 5.2.4 Count percentage difference mark (error)  $\mu_0$  which is obtained through test with reference value
- 5.3 Experiment of Magnetic Field strength in Solenoid with Density Coil Relationship.
  - 5.3.1 Create a graph between the solenoid coil density as the Y axis and the magnetic field strength as the X axis (first convert the units to SI form)
  - 5.3.2 How is the magnetic field in solenoid a related to the density of the solenoid coils ?
  - 5.3.3 Perform linear regression, then determine  $\mu_0$ .
  - 5.3.4 Count percentage difference mark (error)  $\mu_0$  which is obtained through test with reference value
- 5.4 Experiment on the distribution of magnetic field values inside and outside the solenoid.
  - 5.4.1 How is the distribution of the magnetic field inside the solenoid ? Explain, and correlate the explanation with Ampere's law.
  - 5.4.2 How is the magnitude of the magnetic field outside the solenoid? Explain, and correlate the explanation with Ampere's law.

## 5.5 REFERENCES

- Bertrand, Peggy Ann. (2015) : *AP Physics 1 and 2 Induiry-Based Lab Investigations*, 243 – 245.
- Campbell, Larry., Spector, Donald., Allen, Ted. (2003) : *Physics 160 Laboratory Manual*, Hobart and William Smith Colleges, xix – xxi.
- Giri, P. K., (2005) : *Physics Laboratory Manual for Engineering Undergraduates*, Indian Institute of Technology Guwahati, 107 – 109.
- Halliday, D., Resnick, R., Walker, J. (2006) : *Fundamentals of Physics 6<sup>th</sup> Edition*, John Wiley & Sons, 785 – 786.
- Lenda, Andrzej. (2000) : *Physics Laboratory*, International School of Technology of Cracow, 41-1 – 41-3.
- Sutrisno, (2000) : *Seri Fisika Dasar*, Penerbit ITB, 84 – 85.